



HEALTHCARE DISINFECTING WIPES

Disinfect in 2 minutes²



- One Step Disinfectant^{3,5}/Deodorizer/Virucidal*/Bactericidal
- Broad spectrum hospital cleaner disinfectant⁵
- Kills 99.9% of SARS-Co V-2 in 1 minute
- Fast drying, leaves no dull or sticky residue
- Kills cold & flu viruses*



HOSPITAL DISINFECTANT

Cross-contamination between treated hard, nonporous surfaces is of major concern. This product has been formulated to aid in the reduction of cross-contamination on treated surfaces³. Use in hospitals, hospices, long-term care facilities, ambulatory care centers, emergency rooms, operating rooms and more.

2-MINUTE CONTACT TIME²

This product kills bacteria and viruses* in just 2 minutes, Norwalk Virus and Feline Calicivirus in 5 minutes, and SARS-Related Coronavirus 2 (SARS-CoV-2) (cause of COVID-19) in 1 minute. See back page for a complete list of contact times.

DISINFECT HASSLE-FREE

Cleanslate Healthcare Disinfecting Wipes is a one-step hospital disinfectant that is effective against a broad spectrum of bacteria when used as directed^{3,5}. Removes and cleans blood, dried on blood, body oils, fecal matter, urine, organic matter, sticky food spills, splatters and messes.

LIST N DISINFECTANT

EPA Reg. #6836-372-50757
Cleanslate Healthcare Disinfecting Wipes are EPA Registered and on the EPA's List N: Disinfectants for Coronavirus (COVID-19).

THIS PRODUCT KILLS:

Bacteria:

Contact Time (min):

Pseudomonas aeruginosa (Pseudomonas)	2
Salmonella enterica (Salmonella)	2
Staphylococcus aureus (Staph)	2
Acinetobacter baumannii	2
Acinetobacter baumannii, Multi-Drug Resistant	2
Burkholderia cepacia	2
Campylobacter jejuni	2
Klebsiella aerogenes	2
Enterococcus faecalis (Enterococcus)	2
Enterococcus faecalis – Vancomycin Resistant	2
Escherichia coli [E. coli]	2
Escherichia coli – NDM-1 Positive	2
Escherichia coli O157:H7	2
ESBL Escherichia coli	2
Carbapenem Resistant Klebsiella pneumoniae	2
ESBL Klebsiella pneumoniae	2
Klebsiella pneumoniae (Klebsiella)	2
Klebsiella pneumoniae – NDM-1 Positive	2
Listeria monocytogenes	2
Proteus vulgaris	2
Salmonella enterica serovar Typhi (Salmonella)	2
Serratia marcescens [Serratia]	5
Shigella dysenteriae [Shigella]	2
Staphylococcus aureus – Community Associated Methicillin-Resistant (CA-MRSA) (NRS384)	2
Staphylococcus aureus – Community Associated Methicillin-Resistant (CA-MRSA) (NRS123)	2
Staphylococcus aureus – Methicillin-Resistant (MRSA)	2
Staphylococcus aureus – Multi-Drug Resistant	2
Staphylococcus aureus – Vancomycin Intermediate Resistant	2
Staphylococcus aureus – Vancomycin Resistant (VRSA)	2
Streptococcus pneumoniae – Penicillin Resistant	2
Streptococcus pyogenes associated with scarlet fever	2
Vibrio cholerae	2

Viruses:

*Hepatitis B Virus (HBV)	2
*Hepatitis C Virus (HCV)	2
*Herpes Simplex Virus Type 1	2
*Herpes Simplex Virus Type 2	2
*HIV-1 (AIDS Virus)	2
*Human Coronavirus	2
*Human Influenza A Virus (H7N9)	2
*Influenza A Virus H3N2	2
*Influenza A Virus	2
*2009-H1N1 Influenza A Virus (Novel H1N1)	2
*Norwalk Virus (Norovirus)	5
*Respiratory Syncytial Virus (RSV)	2
*Rotavirus (Strain WA)	2
*SARS Associated Coronavirus(SARS)	2
*SARS-Related Coronavirus 2 (SARS-CoV-2) (cause of COVID-19)	1

Animal Viruses:

*Avian Influenza Virus (H3N2)	2
*Avian Influenza Virus (H5N1)	2
*Feline Calicivirus	5

Fungi:

Candida albicans	2
------------------	---

USE IN:

Hospitals, hospices, long-term care facilities, nursing homes, dental offices, medical offices, home health care, laboratories, examination rooms, ambulatory care centers, emergency rooms, operating rooms, anesthesia rooms, CAT labs, ICU areas, critical care units, nurseries, nurse's stations, out-patient facilities, patient room areas, recovery rooms, rehabilitation centers, public restrooms, pharmacies, laundry rooms, locker rooms, bathrooms, institutional facilities, and waiting rooms.



USE ON WASHABLE HARD, NONPOROUS SURFACES OF:

Bedside tables, cabinets, counters, cupboards, chairs, medication carts, bedpans, charging stations, diagnostic equipment, docking stations, equipment tables and surfaces, headboards, footboards, side rails, stools, stretchers, spine backboards, instrument stands, lab benches, nurse call devices, paddles, phlebotomy trays, scales, patient monitoring equipment, car seats, diaper changing tables, diaper pails, highchairs, bathroom fixtures, closet handles, faucets, showers, bathtubs, toilet seats, light switches, doorknobs, floors, and walls.

NON-CRITICAL MEDICAL DEVICE SURFACES OF:

Crutches, defibrillators, gurneys, non-porous mattresses, examining tables, plastic mattress covers, remote controls, rescue tools, resuscitators, stretchers, stethoscopes, walkers, wheel chairs, ultra-sound transducers and equipment.

This product is not for use on medical device surfaces.

² Excluding (Norwalk Virus, Feline Calicivirus, SARS-Related Coronavirus 2 (SARS-CoV-2) (cause of COVID-19).

³ Hard, non-food, non-porous surfaces

⁵ When used according to the directions for disinfection

Carefully read the label and only use the product as described in its directions.